

## dsDNase, Double-strand specific DNases

Cat #: C-BSM0206

Size: 50 rxns

Storage: -20°C

### Components

Components	Size
dsDNase(1 rxn/ $\mu$ l)	50 $\mu$ l
10 $\times$ dsDNase Buffer	200 $\mu$ l

### Product Introduction

dsDNase is a type of endonuclease that cleaves the phosphodiester bonds in DNA, generating oligonucleotides with a 5'-phosphate and a 3'-hydroxyl group. dsDNase specifically digests double-stranded DNA (dsDNA) without digesting single-stranded DNA, primers, probes, and RNA. dsDNase is thermosensitive and rapidly inactivated at 55°C. It is primarily used to quickly remove genomic DNA contamination from RNA samples prior to reverse transcription experiments. Compared to the traditional method of using DNase I to remove genomic DNA contamination, dsDNase does not require additional EDTA for inactivation, reducing RNA damage, saving experimental time, and ensuring accurate quantification of RNA levels.

### Enzyme Activity Definition

According to the Kunitz assay method, under pH 5.0 at 25°C, the enzyme activity is defined as 1 unit (U) when it causes an increase in absorbance of 0.001 at 260 nm per minute using an excess amount of high molecular weight DNA as the substrate.

### Storage Buffer

25 mM Tris-HCl (pH 7.5, at 25°C), 2.0 mM MgCl<sub>2</sub>, 10 mM NaCl, 0.01% (v/v) Triton X-100, 50% (v/v) glycerol.

## Inhibition and Inactivation

**Inhibition:** Metal ions, EDTA, SDS, DTT,  $\beta$ -mercaptoethanol, high salt ion concentrations, etc., can inhibit the activity of dsDNase.

**Inactivation:** Incubation at 55°C for 5 minutes.

## Quality Control

### Protein purity

Determined by SDS-PAGE and Coomassie Brilliant Blue staining, the enzyme purity is  $\geq 90\%$ .

### RNA nuclease activity

The enzyme is incubated with RNA substrate at 37°C for 1 hour, and gel electrophoresis is performed to confirm no degradation of the RNA substrate.

### Functional testing

The product was tested for the removal of genomic DNA from RNA samples and subsequent RT-qPCR amplification, showing a removal rate of  $\geq 99.9\%$ . The quantity of RNA is not affected by dsDNase treatment.

## Instructions for use:

1. Prepare the reaction mixture on ice as follows:

Components	Volume
dsDNase	1 $\mu$ l
10 $\times$ dsDNase Buffer	10 $\mu$ l
Template RNA:	X $\mu$ l
Total RNA	1 pg~5 $\mu$ g/10 $\mu$ l
mRNA	0.1 pg~500 ng/10 $\mu$ l
Specific RNA	0.01 ng~500 ng/10 $\mu$ l
Nuclease-Free Water	To 10 $\mu$ l

2. Gently tap and mix the reaction mixture, then incubate the mixture at 37°C for 2-5 minutes.
3. Heat inactivate at 65°C for 2 minutes, and quickly place the obtained RNA on ice for subsequent experiments.

**For long-term storage, store at -80°C, avoiding repeated freeze-thaw cycles.**

## Notes

1. If the downstream application of the RNA sample is RT-PCR and the target gene length is  $\geq 3$  kb, add a final concentration of 10 mM DTT before the inactivation step.
2. To prevent RNA degradation, add an appropriate amount of RNase inhibitor to the reaction mixture.

## Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.