

Heparin Sodium Anticoagulant Rabbit Blood

Cat #: A-CSH058

Size: 100mL

Storage: 2-8 °C (1 month)

Product Description

Heparin Sodium Anticoagulant Rabbit Blood is prepared by aseptically collecting blood from rabbits and subjecting it to heparin sodium anticoagulation treatment.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Avoid freeze-thaw cycles.
2. Gently invert and mix daily.
3. Handle aseptically, store at low temperature, and use promptly after opening.
4. Avoid prolonged or long-distance transportation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.