

Beagle Plasma (Sodium Citrate Anticoagulant)

Cat #: A-CSH068

Size: 100mL

Storage: -20 °C (1 year)

Product Description

Beagle Plasma (Sodium Citrate Anticoagulant) is prepared by aseptically collecting sodium citrate anticoagulated blood from Beagle dogs and obtaining the supernatant through centrifugation.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Please ensure aseptic procedures and avoid repeated freeze-thaw cycles.
2. The flocculent precipitate that may appear after plasma freeze-thawing is typically caused by the denaturation of lipoproteins during the freezing and thawing process. It does not affect the quality of the product and can be removed by centrifugation at 3500rpm for 10 minutes.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.