

Pig Plasma (Sodium Citrate Anticoagulant)

Cat #: A-CSH108

Size: 100mL

Storage: -20 °C (1 year)

Product Description

Porcine plasma (anticoagulated with sodium citrate) is prepared by aseptically collecting porcine blood anticoagulated with sodium citrate and obtaining the supernatant through centrifugation.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Please ensure aseptic procedures and avoid repeated freeze-thaw cycles.
2. The flocculent precipitate that may appear after plasma freeze-thawing is typically caused by the denaturation of lipoproteins during the freezing and thawing process. It does not affect the quality of the product and can be removed by centrifugation at 3500rpm for 10 minutes.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.