

Normal Calf Serum

Cat #: A-CSH221

Size: 100 mL, 500 mL

Storage: -20°C (5 years)

Product Description

Normal Calf Serum is obtained by sterile collection of calf blood, coagulation and centrifugation.

Notes

1. Use according to specific experimental requirements.
2. Avoid repeated freeze-thaw cycles. Repeated freeze-thaw cycles may result in precipitation, which can be removed by centrifugation at 3500 rpm for 10 minutes under sterile conditions.
3. Exercise caution and maintain sterile conditions during handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.