

Normal Guinea Pig Serum

Cat #: A-CSH229

Size: 100mL, 500mL

Storage: -20 °C (5 years)

Product Description

Guinea pig serum is obtained by aseptically collecting blood from guinea pigs and obtaining the serum after coagulation and centrifugation.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Take care to avoid repeated freeze-thaw cycles. Repeated freeze-thaw cycles can lead to precipitation. If precipitation occurs, it can be removed by centrifuging at 3500 rpm for 10 minutes under aseptic conditions.
2. Pay attention to aseptic techniques during the operation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.