

## Normal Balb/C Mouse Serum

Cat #: A-CSH235

Size: 100mL, 500mL

Storage: -20 °C (5 years)

### Product Description

Babysitter mouse serum is obtained by aseptically collecting blood from Babysitter mice and obtaining the serum after coagulation and centrifugation.

### Usage

It should be used according to the specific experiment requirements.

### Notes

1. Take care to avoid repeated freeze-thaw cycles. Repeated freeze-thaw cycles can lead to precipitation. If precipitation occurs, it can be removed by centrifuging at 3500 rpm for 10 minutes under aseptic conditions.
2. Pay attention to aseptic techniques during the operation.

### Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.