

Sterile Defibrinated Sheep Blood

Cat #: A-CSH253

Size: 100mL, 500mL, 1000mL

Storage: 2-8 °C (1 month)

Product Description

Aseptically defibrinated sheep blood is prepared by aseptically collecting blood from sheep and subjecting it to defibrination.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Avoid freeze-thaw cycles.
2. Gently invert and mix daily.
3. Handle aseptically, store at low temperature, and use promptly after opening.
4. Avoid prolonged and long-distance transportation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.