

1× PBST (0.01M, pH 7.2-7.4)

Cat #: A-CSH487

Size: 500mL

Storage: Room temperature (12 months)

Product Description

PBST, also known as PBS-Tween, is composed of sodium chloride, disodium hydrogen phosphate, and monosodium dihydrogen phosphate. It contains 0.05% Tween 20 and does not contain Ca^{2+} or Mg^{2+} . It is commonly used as a washing solution for ELISA and can also be used to clean PVDF membranes or NC membranes in Western Blot experiments.

Usage

Follow the specific experimental requirements.

Notes

For your safety and health, please wear lab coat and disposable gloves during operation.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.