

## SDS Solution (20%)

Cat #: A-CSH490

Size: 100mL

Storage: Room temperature (12 months)

### Product Description

SDS, also known as Sodium dodecyl sulfate, is a commonly used biochemical reagent with a wide range of applications, including SDS-PAGE and bacterial alkaline lysis.

### Usage

Follow the specific experimental requirements.

### Notes

1. If the ambient temperature is low, precipitation may occur in the solution. You can use a water bath to heat and dissolve it.
2. SDS solution can cause irritation to the human body, please handle with care.
3. For your safety and health, please wear lab coat and disposable gloves during operation.

### Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.