

Glycine Solution (1mol/L)

Cat #: A-CSH491

Size: 100mL

Storage: Room temperature (12 months)

Product Description

Glycine, also known as aminoacetic acid, has a molecular formula of $\text{NH}_2\text{CH}_2\text{COOH}$, a molecular weight of 75.07, and CAS number 56-40-6. The Glycine Solution (1mol/L) is prepared by mixing 1 mole of glycine with distilled water.

Usage

Follow the specific experimental requirements.

Notes

For your safety and health, please wear lab coat and disposable gloves during operation.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.