

Sodium Dihydrogen Phosphate-Citric Acid Buffer (pH 2.2-8.0)

Cat #: A-CSH552

Size: 500mL

Storage: 2~8°C (3 months)

Product Description

Sodium Dihydrogen Phosphate-Citric Acid Buffer (pH 2.2~8.0) is prepared by mixing 0.2M Sodium Dihydrogen Phosphate and 0.1M Citric Acid in different proportions. Customers can choose from pH values ranging from 3.0 to 6.6 according to their specific needs.

Usage

Follow the specific requirements of your experiment.

Notes

For your safety and health, please wear lab attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.