

EDTA Solution (1mol/L, pH 8.0)

Cat #: A-CSH622

Size: 100mL, 500mL

Storage: Room temperature (12 months)

Product Description

EDTA Solution (1M, pH 8.0) is prepared using EDTA•2Na, adjusted to pH 8.0, and subjected to high-pressure sterilization.

This solution is commonly used for chelating metal ions and is a commonly used reagent in molecular biology.

Usage

Operate according to specific experimental requirements.

Notes

1. When sterility is required for the experiment, sterile techniques should be employed.
2. For your safety and health, please wear laboratory attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.