

## Calcium Chloride Solution, $\text{CaCl}_2$ (4.5mol/L, Sterile)

Cat #: A-CSH655

Size: 100mL

Storage: Room temperature (6 months)

### Product Description

Calcium chloride is a fundamental component of various research reagents and can also be used as a desiccant. Calcium chloride solution (4.5 mol/L, sterile) is composed of calcium chloride and deionized water, and it undergoes high-pressure sterilization.

### Usage

Operate according to specific experimental requirements.

### Notes

1. Pay attention to aseptic techniques to avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves during the operation.

### Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.