

Krebs-Ringer-Phosphate-HEPES, KRPH Buffer

Cat #: A-CSH683

Size: 500mL

Storage: 2~8°C (3 months)

Product Description

Krebs-Ringer-Phosphate-HEPES (KRPH) Buffer contains HEPES, potassium dihydrogen phosphate, magnesium sulfate, calcium chloride, sodium chloride, potassium chloride, and others. It has a pH value of 7.4. Similar to Ringer's solution, it is used for in vitro experiments on animals, perfusion of muscles, tissue washing, and maintaining normal physiological functions of isolated tissues.

Usage

Operate according to the specific requirements of the experiment.

Notes

For your safety and health, please wear laboratory clothing and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.