

Potassium Chloride Solution, KCl (3.8mol/L, Sterile)

Cat #: A-CSH702

Size: 100mL

Storage: Room temperature (12 months)

Product Description

Potassium chloride is a fundamental component of various research reagents and can be used to regulate electrolyte balance in the human body. The potassium chloride solution (3.8 mol/L, sterile) is composed of potassium chloride and deionized water, and has undergone high-pressure sterilization.

Usage

Operate according to the specific requirements of the experiment.

Notes

1. Pay attention to aseptic operation to avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.