

DMEM Medium (High Glucose)

Cat #: A-CSH771

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

DMEM (Dulbecco's Modified Eagle Medium) is a widely used basal medium that supports the growth of various mammalian cells, including primary fibroblasts, neurons, glial cells, huvecs, and smooth muscle cells, as well as cell lines such as HeLa, 293, Cos-7, and PC-12. In addition to the incomplete DMEM medium (high sugar), we also provide DMEM medium with added dual antibiotics and fetal bovine serum to meet different customer needs.

This medium contains: 4.5 g/L glucose, 0.584 g/L L-glutamine, phenol red, sodium pyruvate, sodium bicarbonate, etc.

This medium does not contain HEPES.

Notes

1. When using this product for specific experiments, it is important to take precautions to avoid microbial contamination, and avoid freeze-thaw cycles.
2. This product is intended for scientific research purposes by professionals only and is not to be used for clinical diagnosis or treatment, or for food or drug purposes.
3. For your safety and well-being, please wear laboratory attire and disposable gloves when handling this product.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.