

Ham'S F-12K Medium (With Penicillin-Streptomycin)

Cat #: A-CSH773

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

Ham's F-12K (Ham's F-12K (Kaighn's) Medium) is an improved version of Ham's F-12K, which increases the concentration of amino acids and pyruvate, reduces the amount of glucose, and modifies the composition and content of salts based on Ham's F-12. Ham's F-12K medium was originally designed for the culture of primary human liver cells as well as differentiated cells from rats and chickens. In addition to the incomplete Ham's F-12K medium, we also provide Ham's F-12K medium supplemented with antibiotics and fetal bovine serum to meet different customer needs. This medium contains: 1.26 g/L glucose, 0.292g/L L-glutamine, phenol red, sodium bicarbonate, 50KU penicillin, 50mg streptomycin, etc. This medium does not contain HEPES.

Usage

Use according to specific experimental requirements and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination and freeze-thaw cycles.
2. This product is intended for scientific research use by professionals only and is not for clinical diagnosis or treatment, or for use in food or drugs.
3. For your safety and health, please wear a lab coat and disposable gloves during operation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.