

DMEM Medium (No Glucose, Without Phenol Red , With Penicillin-Streptomycin)

Cat #: A-CSH799

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

DMEM (Dulbecco's Modified Eagle Medium) is a widely used basal medium designed to support the growth of various mammalian cells, including primary fibroblasts, neurons, glial cells, HUVECs, smooth muscle cells, as well as cell lines such as HeLa, 293, COS-7, and PC-12.

This medium contains 0.584g/L L-glutamine, sodium pyruvate, sodium bicarbonate, 50KU penicillin, and 50mg streptomycin. It should be protected from microbial contamination and avoid repeated freezing and thawing.

Usage

Use according to specific experimental requirements, and take precautions to avoid microbial contamination.

Notes

1. Take precautions to prevent microbial contamination and avoid repeated freezing and thawing.
2. This product is intended for scientific research purposes by professionals only and should not be used for clinical diagnosis or treatment, or for food or drugs.
3. For your safety and health, wear laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.