

IL27 Monoclonal Antibody

Overview

Catalog No.	BGT-ANT-12770
Host species	Mouse
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	ELISA, FCM
Target	p28, Interleukin-27 subunit alpha, IL-27 subunit alpha, IL27A, IL27, IL-27-A, Interleukin-30, IL30, IL27-A, IL-27B, Interleukin-27 subunit beta, EBI3, Epstein-Barr virus-induced gene 3 protein, EBV-induced gene 3 protein, IL27B, IL-27 subunit beta
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	Q8NEV9 & Q14213
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Clone ID	SAA0393
Note	For research use only.