

PEDV N / Nucleoprotein Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-13709
Alternative Names	Nucleoprotein, Nucleocapsid protein, NC, Protein N, N, 6, Porcine transmissible gastroenteritis coronavirus (strain Purdue) (TGEV)
Host species	Rabbit
Species reactivity	Porcine epidemic diarrhea virus (PEDV)
Immunogen	E. coli - derived recombinant PEDV N/Nucleoprotein (Pro15-Asn441).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, WB
Purification	Antigen affinity purification.
Accession	B7UCM9
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.