

Pseudomonas aeruginosa PcrV Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-13946
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, WB: 1:1000-1:4000
Species reactivity	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Immunogen	E. coli - derived recombinant Pseudomonas aeruginosa PcrV (Leu132-Ile294).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	Type III secretion protein PcrV, pcrV, PA1706
Purification	Purified by antigen affinity column.
Accession	G3XD49
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.