

Newborn Calf Serum

Cat #: A-CSH012

Size: 100mL, 500mL

Storage: -20 °C (5 years)

Product Description

Newborn calf serum is obtained by aseptically collecting blood from newborn calves and obtaining serum after coagulation and centrifugation.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Avoid repeated freeze-thaw cycles, as they may result in precipitation. If precipitation occurs, it can be removed by centrifuging at 3500 rpm for 10 minutes under sterile conditions.
2. Pay attention to aseptic techniques during handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.