

Heparin Sodium Anticoagulant Bovine Blood

Cat #: A-CSH018

Size: 100mL

Storage: 2-8 °C (1 month)

Product Description

Heparin sodium anticoagulated bovine blood is prepared by aseptically collecting bovine blood and treating it with heparin sodium anticoagulant.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Avoid freeze-thaw cycles.
2. Gently invert and mix daily.
3. Handle aseptically, store at low temperature, and use promptly after opening.
4. Avoid prolonged or long-distance transportation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.