

# **EZMeta™ Total Cholesterol (TC) Colorimetric Assay Kit**

Cat #: D-AKC2220

Size: 48T / 96T

Storage: Stored at 4°C for 12 months, protected from light

#### **Product Information**

Applicable samples: Serum, Plasma, Animal Tissues, Cells, Bacteria

# **Assay Principle**

Total Cholesterol (TC) includes free cholesterol and cholesterol esters. Tissue TC refers to the sum of cholesterol contained in all lipoproteins in the tissue. EZMeta™ Total Cholesterol (TC) Colorimetric Assay Kit provides a simple method for detecting TC concentration in a variety of biological samples such as serum, plasma, animal tissues and cells, bacteria. In the assay, esterase catalyzes the hydrolysis of cholesterol esters to produce Free Cholesterol (FC) and Free Fatty Acids (FFA). Then Cholesterol oxidase catalyzes FC to produce Δ4-cholestenone and H<sub>2</sub>O<sub>2</sub>. Further, H<sub>2</sub>O<sub>2</sub>, 4-aminoantipyrine and phenol can be catalyzed by peroxidase to form red quinone compounds which has a characteristic absorption peak at 500 nm, the color depth is proportional to the TC content.

# **Materials Supplied and Storage Conditions**

Kit components	Size		Ctourses sounditions	
	48 T	96 T	Storage conditions	
TC Working Solution	10 mL	20 mL	4°C	
Standard	1 mL	1 mL	4°C, protected from light	





## **Materials Required but Not Supplied**

- ·Microplate reader or visible spectrophotometer capable of measuring absorbance at 500 nm
- ·96-well plate or microglass cuvette, precision pipettes, disposable pipette tips
- ·Ice maker, refrigerated centrifuge, water bath
- ·Anhydrous ethanol
- ·Homogenizer (for tissue samples)

#### **Reagent Preparation**

**TC Working Solution:** Ready to use as supplied. Equilibrate to 37°C before use. Store at 4°C.

Standard: Ready to use as supplied. Equilibrate to room temperature before use. Store at 4°C, protected from light.

Standard Curve Setting: Dilute the 2  $\mu$ mol/mL standard with absolute ethanol to 1, 0.5, 0.25, 0.125, 0.063, 0.031  $\mu$ moL/mL as indicated in the table below.

Num.	Volume of Standard	Volume of Anhydrous	The Concentration of Standard	
	volume of Standard	Ethanol (μL)	(μmol/mL)	
Std.1	200 μL of Standard	0	2	
Std.2	100 μL of Std.1 (2 μmol/mL)	100	1	
Std.3	100 μL of Std.2 (1 μmol/mL)	100	0.5	
Std.4	100 μL of Std.3 (0.5 μmol/mL)	100	0.25	
Std.5	100 μL of Std.4 (0.25 μmol/mL)	100	0.125	
Std.6	100 μL of Std.5 (0.125 μmol/mL)	100	0.063	
Std.7	100 μL of Std.6 (0.063 μmol/mL)	100	0.031	

#### **Sample Preparation**

Note: Fresh samples are recommended, If not assayed immediately, samples can be stored at -80°C for one month.

- 1. Animal tissues: Weigh about 0.1 g tissues and add 1 mL anhydrous ethanol. Homogenize on ice. Centrifuge at 8,000 g for 10 min at 4°C. Use supernatant for assay, and place it on ice to be tested.
- 2. Cells or Bacteria: Collect 5×10<sup>6</sup> cells or bacteria into the centrifuge tube, discard the supernatant after centrifugation;





add 1 mL anhydrous ethanol to ultrasonically disrupt the cells or bacteria in ice bath 5 min (power 20% or 200 W, ultrasonic 3 s, interval 7 s, repeat 30 times). Centrifuge at 8,000 g for 10 min at 4°C. Use supernatant for assay, and place it on ice to be tested.

3. Serum (plasma) samples: Directly test.

Note: It will be better to quantify the total protein with EZMeta™ Total Cholesterol (TC) Colorimetric Assay Kit, if the content is calculated by protein concentration.

# **Assay Procedure**

- 1. Preheat the microplate reader or visible spectrophotometer for more than 30 min, and adjust the wavelength to 500 nm, visible spectrophotometer was returned to zero with deionized water.
- 2. Warm TC Working Solution to 37°C for more than 30 min before use.
- 3. Add the following reagents to the 96-well plate or microglass cuvette:

Reagent	Blank Well (μL)	Standard Well (μL)	Test Well (μL)
Anhydrous Ethanol	50	0	0
Standard	0	50	0
Sample	0	0	50
TC Working Solution	150	150	150

4. Mix well and then incubatefor 15 min at room temperature. Then reading the values at 500 nm. The absorbance of blank well, standard well, test well recorded as  $A_{Blank}$ ,  $A_{Standard}$  and  $A_{Test}$ . Finally, calculate  $\Delta A_{Test} = A_{Test} - A_{Blank}$ ,  $\Delta A_{Standard} = A_{Standard} - A_{Blank}$ .

Note: Blank well and standard well only need to measure 1 time. In order to guarantee the accuracy of experimental results, need to do a pre-experiment with 2-3 samples. If  $\Delta A_{Test}$  is less than 0.001, increase the sample quantity appropriately. If  $\Delta A_{Test}$  is greater than 2.0, the sample can be appropriately diluted with anhydrous ethanol, the calculated result multiplied by the dilution factor, or decrease the sample quantity appropriately. The test time should not exceed 1 h.





## **Data Analysis**

Note: We provide you with calculation formulae, including the derivation process and final formula. The two are exactly equal. It is suggested that the concise calculation formula in **bold** is final formula.

1. Drawing of standard curve

With the concentration of the Standard solution as the y-axis and the  $\Delta A_{Standard}$  as the x-axis, draw the standard curve. Substitute the  $\Delta A_{Test}$  into the equation to obtain the y value ( $\mu moL/mL$ ).

- 2. Calculate the content of TC in sample
- (1) By sample fresh weight

TC ( $\mu$ mol/g)=y÷(W÷V<sub>Extraction</sub>)×n=y÷W×n

(2) Calculated by protein concentration

TC (μmol/mg prot)=y÷Cpr×n

(3) Calculated by cells or bacteria number

TC (μmol/10<sup>4</sup>)=y÷(cells or bacteria number÷V<sub>Extraction</sub>)×n=y÷500×n

(4) Calculated by liquid volume

TC ( $\mu$ mol/mL)= $y \times n$ 

Where: W: sample weight, g; V<sub>Extraction</sub>: Anhydrous ethanol volume added, 1 mL; n: dilution factor; Cpr: supernatant protein concentration, mg/mL; 500: Total number of cells or bacteria, 5×10<sup>6</sup>.

# **Typical Data**

Typical standard curve:








Figure 1. Standard curve for TC.

#### Examples




Figure 2. TC content in animal tissues. Assays were performed following kit protocol.

#### **Disclaimer**

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.

