

IHCAb™ Calcitonin mouse mAb (BGT124)

Cat #: B-IMW6824

Size: 100 µL

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Product Information

Applications/Dilution: IHC-p 1:100-500, IF 1:100-500

Isotype/Source: Mouse, Monoclonal/IgG1, Kappa

Specificity: The antibody can specifically recognize human Calcitonin protein.

Subcellular Location: Secreted

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.21% sodium azide

Storage: Store at -15°C to -25°C




Fig.1. Human thyroid atypical adenoma tissue was stained with anti-Calcitonin antibody




Fig.2. Human thyroid medullary carcinoma tissue was stained with anti-Calcitonin antibody

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.