

IHCAb™ Calcitonin mouse mAb (BGT285)

Cat #: B-IMW6825

Size: 100 µL

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Product Information

Applications/Dilution: IHC-p 1:100-500, IF 1:100-500

Isotype/Source: Mouse, Monoclonal/IgG2a, Kappa

Specificity: The antibody can specifically recognize human Calcitonin protein.

Subcellular Location: Secreted

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.22% sodium azide

Storage: Store at -15°C to -25°C




Fig. Human thyroid atypical adenoma tissue was stained with anti-Calcitonin antibody.

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.