

IHCAb™ Calponin-1 mouse mAb (BGT133)

Cat #: B-IMW6827

Size: 100 µL

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Product Information

Applications/Dilution: IHC-p 1:100-500, WB 1:200-1000, IF 1:100-500

Isotype/Source: Mouse, Monoclonal/IgG2a, Kappa

Specificity: The antibody can specifically recognize human and mouse calponin-1 protein. In western blotting of HeLa, LnCap and mouse muscle lysates, the antibody can label a 34KDa band corresponding to calponin-

Subcellular Location: cytoskeleton, focal adhesion

Expression: Smooth muscle, and tissues containing significant amounts of smooth muscle

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.146% sodium azide

Storage: Store at -15°C to -25°C



Fig. Human appendix tissue was stained with anti-calponin-1 antibody.

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.