

IHCAb™ Lysozyme mouse mAb (BGT191)

Cat #: B-IMW6885

Size: 100 µL

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Product Information

Applications/Dilution: IHC-p 1:100-500, WB 1:200-1000, IF 1:100-500

Isotype/Source: Mouse, Monoclonal/IgG3, Kappa

Specificity: The antibody can specifically recognize human Lysozyme protein

Subcellular Location: Secreted

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.156% sodium azide

Storage: Store at -15°C to -25°C




Fig.1. Human lung tissue was stained with Anti-Lysozyme Antibody




Fig.2. Human malignant fibrous histiocytoma tissue was stained with Anti-Lysozyme Antibody

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.