

## IHCAb™ Mammaglobin mouse mAb (BGT079)

Cat #: B-IMW6888

Size: 100 µL

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

### Product Information

**Applications/Dilution:** IHC-p 1:100-500, IF 1:100-500

**Isotype/Source:** Mouse, Monoclonal/IgG1, Kappa

**Specificity:** The antibody can specifically recognize human Mammaglobin protein

**Subcellular Location:** Secreted

**Expression:** Mammary gland specific. Over-expressed in breast cancer

**Formulation:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.167% sodium azide

**Storage:** Store at -15°C to -25°C




Fig. Human breast tissue was stained with Anti-Mammaglobin Antibody.

### Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.