

IHCAb™ hCG β (BGT184) Mouse mAb

Cat #: B-IMW6959

Size: 100 μ L

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Product Information

Applications/Dilution: IHC-p 1:200-400

Isotype/Source: Mouse, Monoclonal/IgG1, Kappa

Specificity: The antibody can specifically recognize β subunits of human CG, and has no cross reaction with LH, FSH or TSH proteins

Subcellular Location: Cytoplasmic

Expression: Placenta

Formulation: PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein

Storage: Store at -15°C to -25°C




Fig. Human placenta tissue was stained with Anti-hCG β Antibody

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.