

IHCab™ GCDFP-15 (BGT175) Mouse mAb

Cat #: B-IMW6990

Size: 100 µL

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Product Information

Applications/Dilution: IHC-p 1:200-400

Isotype/Source: Mouse, Monoclonal/IgG1, Kappa

Specificity: The antibody can specifically recognize human GCDFP-15 protein

Subcellular Location: Cytoplasmic

Expression: Breast

Formulation: PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein

Storage: Store at -15°C to -25°C




Fig. Human breast tissue was stained with Anti-GCDFP-15 Antibody

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.