

**Bovine Serum Albumin(IgG, Protease, Fatty-acid free),
Premium Grade**

Catalog Number	BGT-CMP-104
Level	Premium Grade, Fatty-acid free
Sizes	100g, 1000g, 10kg
Appearance	Lyophilized, White to light yellow powder
Source	Normal, healthy bovine
Purification	Heat-Shock
Moisture	< 2.0%
Ash	< 1.0%
pH	6.8-7.2 (10% solution in H ₂ O)
Purity	≥ 99% (SDS-PAGE)
Protein	≥ 99% (Biuret method)
Solubility	4% solution in H ₂ O, dissolution time ≤15mins
Protease	None Detected
Endotoxin	≤ 1EU/mg
Bioburden	≤ 10cfu/g
BlgG	None Detected
NEFA	< 0.01%
Storage	2-8℃