

GMP Grade, Recombinant Human Interleukin-2 (rHuIL-2)

Cat#: C-BETA001

Size: 100µg, 1mg

Source: *Escherichia coli*.

Product Information

Molecular Weight: Approximately 15.5 kDa, a single non-glycosylated polypeptide chain containing 134 amino acid.

Description : Accession # P60658, Ala21-Thr153 (C125A), with an N-terminal Ser.

SDS-PAGE: 15.5 kDa, reducing conditions

Purity: > 95 %, as determined by SDS-PAGE, under reducing non-reducing conditions, visualized by coomassie staining.

Endotoxin: Less than 0.01 EU/µg of IL-2 as determined by kinetic Limulus Amoebocyte Lysate (LAL) assay.

Biological Activity: Recombinant human IL-2 bioactivity is measured in a cell proliferation assay using CTLL-2 mouse cytotoxic T cells, the EC50 for this effect is 0.183 to 0.243 ng/mL.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation: Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM PB, pH 3.5 with 10% trehalose anhydrous.

Reconstitution: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute to a concentration of 0.1-1.0 mg/mL in **sterile distilled H₂O**. Stock solutions should be apportioned into working aliquots and stored at -20 °C to -70 °C. Further dilutions should be made in appropriate buffered solutions. **Do not reconstitute in cell culture media directly.**


Shipping: 2°C to 8°C. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. A minimum of 12 months from date of shipping when stored at -20°C to -70°C as supplied. 4 weeks at 2°C to 8°C under sterile conditions after reconstitution. 4 months at -20°C to -70°C under sterile conditions after reconstitution.


Usage: Biogradetech IL-2 product can be used for a variety of ex vivo cell culture applications such as research or further manufacturing.

Quality statement: No animal- or human-derived materials were used for the manufacture of this product, unless otherwise stated in the respective Certificate of Origin.


DATA




SDS-PAGE: Recombinant Human IL-2 Protein SDS-PAGE 1 µg/lane of Recombinant Human IL-2(Catalog# C-BETA001) was resolved with SDS-PAGE under reducing (R) and non-reducing (NR) conditions visualized by coomassie staining showing a single band at 15 kDa.



Bioactivity-ELISA: Immobilized Recombinant Human IL-2(Catalog# C-BETA001) at 0.2 µg/well can bind Human IL-2 R alpha, Fc Tag, with a linear range of 378.4 to 508.8 ng/mL.



Bioactivity-BLI: Loaded Human IL-2 R alpha, can bind Recombinant Human IL-2(Catalog# C-BETA001) with an affinity constant of 4.55 nM as determined in BLI assay (Octet Red96e) .



Bioactivity-Cell based assay: Recombinant Human IL-2(Catalog# C-BETA001)stimulates proliferation of CTLL-2 cells,the EC50 for this effect is 0.183 to 0.243 ng/mL.

References

1. Liao W, et al. 2013. Immunity. 38:13-25.
2. Sakaguchi S, et al. 2008. Cell. 133:775-87.
3. Kryczek, I. et al. 2007. J. Immunol. 178:6730.
4. Cote-Sierra J, et al. 2004. Proc Natl Acad Sci USA. 101:3880-5.

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.