

Anti-DDDDK Tag Rabbit Pab Antibody

Cat #: D-AKE2011

Size: 50 μ L / 200 μ L

Storage: Store at -20°C . Avoid repeated freeze / thaw cycles.

Background

The DYKDDDDK peptide (Flag-tag) is a polypeptide protein tag that can be added to a protein using recombinant DNA technology. It can be used for affinity chromatography, and then used to separate recombinant, over expressed protein from wild-type protein expressed by the host organism. It can also be used in the isolation of protein complexes with multiple subunits.

Product Information

Applications: WB

Suggested starting dilutions are as follows: WB (1:5000).

Isotype: Rabbit IgG

Reactivity: Mammals, Bacteria

Formulation: Liquid

Concentration: 1 mg/mL

Storage: Store at -20°C . Avoid repeated freeze / thaw cycles.

Storage Buffer: Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.




Fig. Western blot analysis of 1 μ g Flag fusion protein with Anti-Flag rabbit polyclonal antibody in 1:3000 (lane A) and 1:6000 (lane B) dilutions.

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.