

Biotin, Goat Anti-Rabbit IgG (H+L)

Cat #: D-AKE21220

Size: 100µl / 500µl

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Background

Biogradetech secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

Product Information

Applications: WB, IHC, ELISA

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:1000-1:10000), IHC (1:500-1:1000), ELISA (1:1000-1:10000).

Isotype: Goat IgG

Reactivity: Rabbit

Formulation: Liquid

Concentration: 1 mg/ml

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Storage Buffer: Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.