

IF-Light™ Green Donkey Anti-Rabbit IgG (H+L)

Cat #: D-AKE24221

Size: 100µl / 500µl

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Background

Biogradetech secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application.

Product Information

Applications: IF, ICC, FCM

Suggested starting 1:50-1:1000 dilutions for most fluorescent applications.

Isotype: Donkey IgG

Reactivity: Rabbit

Formulation: Liquid

Concentration: 0.5 mg/ml

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

Storage Buffer: Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.