

## IF-Light™ Orange Donkey Anti-Mouse IgG (H+L)

Cat #: D-AKE24311

Size: 100µl / 500µl

Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.

### Background

Biogradetech secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application.

### Product Information

**Applications:** IF, ICC, FCM

Suggested starting 1:50-1:1000 dilutions for most fluorescent applications.

**Isotype:** Donkey IgG

**Reactivity:** Mouse

**Formulation:** Liquid

**Concentration:** 0.5 mg/ml

**Storage:** Store at -20°C. Avoid repeated freeze / thaw cycles.

**Storage Buffer:** Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol.

**Storage Instructions:** Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

### Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.