

## Potassium Oxalate Anticoagulant Horse Blood

Cat #: A-CSH046

Size: 100mL

Storage: 2-8 °C (1 month)

### Product Description

Potassium Oxalate Anticoagulant Horse Blood is prepared by aseptically collecting horse blood and treating it with potassium oxalate anticoagulant.

### Usage

It should be used according to the specific experiment requirements.

### Notes

1. Avoid freeze-thaw cycles.
2. Gently invert and mix daily.
3. Handle aseptically, store at low temperature, and use promptly after opening.
4. Avoid prolonged or long-distance transportation.

### Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.