

Potassium Oxalate Anticoagulant Porcine Blood

Cat #: A-CSH055

Size: 100mL

Storage: 2-8 °C (1 month)

Product Description

Sodium Citrate Anticoagulant Porcine Blood is prepared by aseptically collecting blood from pigs and subjecting it to potassium oxalate anticoagulation treatment.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Avoid freeze-thaw cycles.
2. Gently invert and mix daily.
3. Handle aseptically, store at low temperature, and use promptly after opening.
4. Avoid prolonged or long-distance transportation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.