

## Heparin Sodium Anticoagulant Sheep Blood

Cat #: A-CSH062

Size: 100mL, 500mL, 1000mL

Storage: 2-8 °C (1 month)

### Product Description

Heparin Sodium Anticoagulant Sheep Blood is prepared by aseptically collecting blood from sheep and subjecting it to heparin sodium anticoagulation treatment.

### Usage

It should be used according to the specific experiment requirements.

### Notes

1. Avoid freeze-thaw cycles.
2. Gently invert and mix daily.
3. Handle aseptically, store at low temperature, and use promptly after opening.
4. Avoid prolonged or long-distance transportation.

### Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.