

## Cynomolgus Monkey Plasma (Heparin Sodium Anticoagulant)

Cat #: A-CSH081

Size: 50 mL

Storage: -20°C (12 months)

### Product Description

Cynomolgus Monkey Plasma (Heparin Sodium Anticoagulant) is prepared by aseptically collecting heparin sodium anticoagulated blood from cynomolgus monkeys, followed by centrifugation to obtain the supernatant.

### Notes

1. Use according to specific experimental requirements.
2. Please ensure sterile conditions during handling and avoid repeated freeze-thaw cycles.
3. The flocculent precipitate that may appear after plasma thawing is typically caused by the denaturation of lipoproteins during the freeze-thaw process. It does not affect the quality of the product and can be removed by centrifugation at 3500 rpm for 10 minutes.

### Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.