

Cynomolgus Monkey Plasma (Sodium Citrate Anticoagulant)

Cat #: A-CSH083

Size: 50 mL, 100 mL

Storage: -20°C (12 months)

Product Description

Cynomolgus Monkey Plasma (Sodium Citrate Anticoagulant) is prepared by aseptically collecting sodium citrate anticoagulated blood from cynomolgus monkeys, followed by centrifugation to obtain the supernatant.

Notes

1. Use according to specific experimental requirements.
2. Please ensure sterile conditions during handling and avoid repeated freeze-thaw cycles.
3. The flocculent precipitate that may appear after plasma thawing is typically caused by the denaturation of lipoproteins during the freeze-thaw process. It does not affect the quality of the product and can be removed by centrifugation at 3500 rpm for 10 minutes.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.