

Mouse Plasma PRP (Sodium Citrate Anticoagulant)

Cat #: A-CSH103

Size: 50mL

Storage: 2~8°C (1 month)

Product Description

Mouse plasma PRP (sodium citrate anticoagulant) was prepared by aseptically collecting sodium citrate anticoagulant mouse blood, centrifuging at low speed and taking supernatant.

Usage

It should be used according to the specific experiment requirements.

Notes

Please ensure aseptic procedures and avoid repeated freeze-thaw cycles.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.