

10% Aldehyde-Fixed/Formaldehyde Chicken Red Blood Cells

Cat #: A-CSH140

Size: 50ml 100mL

Storage: 2-8 °C (6 months)

Product Description

Cynomolgus monkey red blood cells are obtained by aseptically collecting anticoagulated blood from cynomolgus monkeys, followed by centrifugation to obtain a red blood cell pellet. The pellet is then washed 2-3 times with an equal volume of PBS until the supernatant becomes clear. Finally, the red blood cells are resuspended in Alsever's solution to a final concentration of 1% cynomolgus monkey red blood cells.

Aldehyde-Fixed/Formaldehyde chicken red blood cells were collected by sterile anticoagulant chicken blood, centrifuged to obtain red blood cell precipitation, washed with 10x PBS 4~5 times, after glutaraldehyde treatment, washed with 10x normal saline 4~5 times. Finally, the red blood cells are resuspended in Alsever's solution to a final concentration of 10% Aldehyde-Fixed/Formaldehyde Chicken Red Blood Cells.

Usage

It should be used according to the specific experiment requirements.

Notes

1. Ensure aseptic handling and store at low temperature. Once opened, please use the product as soon as possible.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.