

Borate Buffer Solution (BBS, pH 9.4)

Cat #: A-CSH428

Size: 500mL

Storage: Room temperature (12 months)

Product Description

Borate Buffer Solution (BBS, pH 9.4) is composed of boric acid, sodium tetraborate, and sodium chloride, with the pH adjusted to 9.4.

Usage

Operate according to the specific requirements of the experiment.

Notes

1. If the amount used each time is small, it can be aliquoted appropriately before use.
2. For your safety and health, wear laboratory attire and use disposable gloves when handling this reagent.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.