

PH Electrode Storage Solution

Cat #: A-CSH444

Size: 100mL

Storage: Room temperature (36 months)

Product Description

Prolonged exposure to acidic or alkaline conditions or improper storage solutions can lead to aging of pH meter electrodes. The pH Meter Electrode Storage Solution effectively protects pH meter electrodes.

Usage

1. Operate according to specific experimental requirements.
2. Generally, thoroughly clean the electrode with ultrapure water and then immerse it in the pH Meter Electrode Storage Solution. Before using the electrode for pH measurement, rinse it thoroughly with ultrapure water.

Notes

For your safety and health, please wear a lab coat and disposable gloves when handling

Disclaimer

For your safety and health, please wear lab attire and disposable gloves while handling.