

Bromophenol Blue Aqueous Solution (1%)

Cat #: A-CSH445

Size: 10mL

Storage: Room temperature (12 months)

Product Description

Bromophenol Blue is a pH indicator that changes color from yellow to blue in the pH range of 3.0 to 4.6. It is commonly used as an electrophoresis tracking dye in gel electrophoresis to monitor the migration of small molecules of nucleic acids or proteins. It is a widely used molecular biology reagent.

Usage

Operate according to specific experimental requirements.

Notes

For your safety and health, please wear a lab coat and disposable gloves when handling

Disclaimer

For your safety and health, please wear lab attire and disposable gloves while handling.