

Iron Chloride Solution (3mol/L)

Cat #: A-CSH449

Size: 100mL

Storage: Room temperature (12 months)

Product Description

Ferric chloride, also known as iron(III) chloride, has a chemical formula FeCl_3 , a molecular weight of 162.2, and a CAS number of 7705-08-0. It can be used for metal etching, wastewater treatment, and as an oxidizing agent.

The Ferric Chloride Solution (3mol/L) is composed of anhydrous ferric chloride and deionized water. This reagent is intended for research purposes only and is not suitable for clinical diagnosis or other applications.

Usage

Follow specific experimental requirements.

Notes

For your safety and health, please wear lab attire and disposable gloves while handling.

Disclaimer

For your safety and health, please wear lab attire and disposable gloves while handling.