

SDS Solution (10%)

Cat #: A-CSH486

Size: 100mL,500mL

Storage: Room temperature (12 months)

Product Description

SDS, also known as sodium dodecyl sulfate, is a commonly used biochemical reagent with a wide range of applications.

It can be used for SDS-PAGE, bacterial lysis, and more.

Usage

Follow the specific experimental requirements.

Notes

1. If the amount used each time is small, it can be appropriately aliquoted before use.
2. For your safety and health, please wear lab coat and disposable gloves during operation.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.